@operSurgical

Endosee[®] System

With Cannula


Sistema Endosee[®]

con cánula

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Avec canule

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@operSurgical

About This Manual

This manual provides instructions on how to use the Endosee® System with the Disposable Cannula safely and effectively for its intended use. It is important to read and observe the information provided in this manual prior to use for proper performance, correct operation, and to ensure patient and operator safety. The user is responsible for operating this device as indicated in this guide.

Endosee System Part Number:

ES9000 (USA, Canada, Mexico):

Contains Display Module, Docking Station, Power Supply (IEC Plug Type A), USB Cable and IFU

Spare Unit Part Numbers*:

ESDISP-SP: Display Module, Spare Unit

ESDOCK-SP: Docking Station, Spare Unit

ESCABL-SP: USB/Power Supply, Spare Units

* Must have previously purchased ES9000 System to order spare units

Endosee Cannula Part Numbers:

ESPX5 - Cannula with endometrial suction pipet curette

ESPY5 - Cannula only

95 Corporate Drive • Trumbull, CT, 06611 • USA

Phone: (800) 243-2974 • Fax: (800) 262-0105

www.coopersurgical.com

Table of Contents

Warranty	5
1 General Description and Intended Use	6
1.1 Indications for Use	7
1.2 Contraindications	7
2 Warnings and Precautions	8
2.1 General.....	8
2.2 Inspection, Use and Disposal.....	9
2.3 Battery and Power Supply.....	10
2.4 Environmental	11
3 Glossary of Symbols	13
4 Preparation for Use	15
4.1 Unpacking and Inspection of Cannula	15
4.2 Power Requirements	15
4.3 Environmental Requirements	16
5 Description of Components	17
5.1 Product Description	17
5.2 Reusable Handheld Display Module	17
5.3 Docking Station	18
5.4 Sterile, Single-use Disposable Diagnostic Cannula.....	19
5.5 Device Classification, Technical and Safety Specification.....	20
5.5.1 Technical Specifications.....	20
6 Basics of Operation and Procedure	21
6.1 Cautions for Handling and Use	21
6.2 Fluid Delivery System.....	21
6.3 Cannula Attachment and removal to the Display Module.....	22
6.4 Power Button Function	23
6.5 Video and Picture (Single Frame Button).....	23
6.6 Touch Screen Functions.....	24
6.7 Patient Examination Procedure.....	30
6.8 Removal of the System	31
6.9 Disposal of Biohazardous Material.....	31
6.10 Battery Charging	31

7	External Interface Information	32
8	EMC Compliance Information	36
8.1	Guidance and Manufacturer's Declaration – Electromagnetics emissions	36
8.2	Guidance and Manufacturer's Declaration – Electromagnetics immunity	37
9	Maintenance	39
9.1	Cleaning/Disinfecting Instructions	39
9.2	Intermediate Level Combined Cleaning/Disinfecting Instructions	39
9.3	Disposable Cannula Handling/Care	41
9.4	Maintenance	41
10	Storage & Shipping	41
11	Technical Assistance	42
12	Troubleshooting Guide	42

Warranty

CooperSurgical, Inc. warrants that the device(s) supplied is free from defects in materials or workmanship.


This warranty is valid only if the product is supplied to the end user by a CooperSurgical, Inc. approved agent or distributor and has been maintained in accordance with procedures documented in the User Manual. If failure occurs from manufacturing defects within 12 months of purchase, CooperSurgical, Inc. will repair or replace the defective item.

1 General Description and Intended Use

Endosee® System with Cannula is intended for use in viewing the cervical canal, uterine cavity, or female urinary tract, including the bladder, for the purpose of performing diagnostic procedures in an outpatient or in an office setting. Hysteroscopy can be used for the evaluation of the uterine cavity for pathology such as submucous myomas, endometrial polyps, malignancy and other conditions.

The Endosee System is a lightweight handheld battery operated portable system. It is used with a single-use Disposable Cannula with a working channel, camera and light source at the distal end to illuminate the area for visualization and image and video capture. The video signal is electronically transferred to the main body of the hysteroscope via an electrical connector. A LCD touch screen Display Module is mounted on the Cannula for viewing.

Note: *The Endosee System cannot be used with an external monitor.*

1. Sterile, Single-use Disposable Cannula
2. Reusable Display Module

Figure 1: Ready-to-use assembly

1.1 Indications for Use

The Endosee® System is used to permit viewing of the adult cervical canal and uterine cavity for the purpose of performing diagnostic procedures. Generally recognized indications for diagnostic hysteroscopy include:

- Abnormal bleeding
- Infertility and pregnancy wastage
- Evaluation of abnormal hysterosalpingogram
- Intrauterine foreign body
- Amenorrhea
- Pelvic pain

The Endosee® System can also be used to permit viewing of the adult urinary bladder through a minimally invasive approach by utilizing natural orifices to access the diagnostic site.

1.2 Contraindications

Hysteroscopy is contraindicated in patients with known or suspected acute pelvic inflammatory disease. Hysteroscopy may be contraindicated by the following conditions depending on their severity or extent:

- Inability to distend the uterus
- Cervical stenosis
- Cervical/Vaginal infection
- Uterine bleeding or menses
- Known pregnancy
- Invasive carcinoma of the cervix
- Recent uterine perforation
- Medical contraindication or intolerance to anesthesia

Cystoscopy is contraindicated in febrile patients with urinary tract infection (UTIs) and those with severe coagulopathy.

2 Warnings and Precautions

This section contains important safety information related to use of the The Endosee® System and the Cannula . Other important safety information is provided throughout this manual in sections that relate specifically to the precautionary information. Read all text surrounding all warning and precautionary information prior to performing any procedures with this equipment.

2.1 General

- **CAUTION:** Federal Law (USA) restricts the device to sale by or on the order of a physician.
- The Endosee® System is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.
- The Endosee® System is for use only by physicians with adequate training in hysteroscopy and/or cystoscopy.
- When a liquid distension medium is used, strict fluid intake and output surveillance should be maintained. Intrauterine instillation exceeding 1 liter should be followed with great care to reduce the possibility of fluid overload.
- Potential complications of continuous flow hysteroscopy:
 - Hyponatremia
 - Hypothermia
 - Uterine perforation resulting in possible injury to bowel, bladder, major blood vessels and ureter
 - Pulmonary edema
 - Cerebral edema
 - Infection
 - Bleeding
 - Pain

The use of normal saline as a distending medium, and limiting the infused volume to less than 1000 mL is recommended to decrease the risk of the above complications.

- Intrauterine distension can usually be accomplished with pressures in the range of 35 - 75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.
- Vaginal ultrasonography before hysteroscopy may identify clinical conditions that will alter patient management.
- Since most urological endoscopy procedures utilize irrigation fluids, it is important that the correct fluid be chosen and that the volumes of inflow and outflow be strictly observed. Appropriate measures should be undertaken to prevent the possibility of fluid overload of the patient

2.2 Inspection, Use and Disposal

- Inspect the integrity of the unit and condition before powering-on the fully assembled Endosee® System. Do not use the system if indications of external damage are observed.
- The Endosee® System contains no operator serviceable components within the enclosure. To avoid electrical shock or damage, the unit must not be disassembled. Doing so will void the warranty.
- Do not operate the Endosee® System if any shipping damage or other defects to the equipment are noted during inspection. Immediately notify CooperSurgical, Inc. Customer Service if any defect is found.
- Do not use the Cannula if the package has been damaged or the expiration date on the label has passed. The cannula is single use only. Do not reuse or resterilize the cannula.
- The Cannula must be disposed of as biohazardous waste according to the safety guidelines of user facility/institution.
- Follow all cleaning procedures provided for the device (applies to reusable Display Module).
- The Endosee® System can only be used with the LCD touch screen display on the Display Module. A separate video monitor cannot be used.

- Do not allow any fluids or substances to be spilled into or around the removable cannula connector on the Display Module, as damage to the electronics may occur.
- Do not misalign the connectors between the removable cannula and the Display Module during assembly, as damage may occur.
- When detaching the cannula from the Display Module, be sure that fluid does not contact the connector in the Display Module. Do not rotate or twist the cannula while removing.
- Do not alter the shape of the distal end of the cannula.
- No portion of the Endosee® System should be opened for cleaning or disinfecting.
- Do not pour water or any cleaning solution directly onto Display Module during cleaning or disinfection. Doing so may damage electronic components and reduce product life.
- Do not immerse any component of the Display Module into solutions during the cleaning or disinfecting processes, as damage to the electronic components may occur.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories other than those provided by CooperSurgical, Inc. could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

2.3 Battery and Power Supply

Comply with the following to prevent damage to the battery:

- Use only the power supply provided with the device. Use of other chargers may damage the device and will void the warranty.

- Charge the battery at room temperature.
- Do not charge the battery near a heat source.
- Do not disassemble or alter the battery in any way.
- Battery is not user-replaceable.
- The Display Module should be fully charged before use. A fully charged battery will provide at least 1 hour of continuous operating time under normal operating conditions.
- Fully charge the Display Module after each procedure. A Low Emission Diode (LED) indicator located on the Docking Station will display a steady green light when the unit is fully charged.
- Refer to battery power indicator on the touch screen Display Module to determine remaining battery power.
- Do not attempt to use any power source or battery not directly provided by CooperSurgical, Inc.

Comply with the following to reduce the risk of burns, fire, electric shock, or injury to persons:

- Do not operate the power supply if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or if it has fallen into water.
- Do not charge the battery where oxygen is being administered.
- Do not unplug power supply by pulling on the cord. To unplug, pull on the plug, not the cord.

2.4 Environmental

- Failure to ensure that all of these environmental requirements are met may result in improper performance or damage to the electronic components within the Endosee® System.
- Due to the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare environments (for example cellular phones, MRI, or active HF surgical equipment, etc.) it is possible that high levels of such interference due to close proximity or strength of a source may result in a disruption of the performance of this device.

- The Endosee® System is not designed for use in environments in which strong interference may impact the performance of the equipment.
- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment
- If this occurs, the site of use should be surveyed to determine the source of this disruption, and the following actions may be taken to eliminate the source:
 - Remove, reorient or relocate the interfering equipment
 - Increase the separation between the interfering equipment and the Endosee® System
 - Incrementally turn off equipment in the vicinity to identify the interfering device
- Do not attempt to transport or ship the Endosee® System without using proper packaging to protect the product.
- This product is intended for indoor use only.

3 Glossary of Symbols

The following symbols are used in this manual, and/or on the product labeling.

The symbols below are referenced in ISO-15223-1:2016

REF. #	Symbol	Title	Description
5.1.1		Manufacturer	Indicates the medical device manufacturer, as defined in EU Directiver 90/385/EECV, 93/42/EEC and 98/79/EC.
5.1.4		Use-by Date	Indicates the date after which the medical device is not to be used.
5.1.5		Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
5.1.6		Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
5.1.7		Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
5.2.3		Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
5.2.4		Do not re-use	Indicates a medical device that is intended for one use, or use on a single patient during a single procedure.
5.2.6		Do not re-sterilize	Indicates a medical device that is not to re-sterilize
5.2.7		Non-sterile	Indicates a medical devicethat has not been subjected to a sterilization process.

REF. #	Symbol	Title	Description
5.2.8		Do not use if package is damaged	Indicates a medical device that should not be used if the package has been broken or damaged.
5.4.4		Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Symbols not listed in the ISO-15223-1:

MEDICAL EQUIPMENT WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARD ONLY IN ACCORDANCE WITH UL60601-1 AND CAN/CSA C22.2 No.601.1

4 Preparation for Use

4.1 Unpacking and Inspection of the Cannula

Follow below to unpack, inspect, and install the Cannula and the Display Module. Each single- use, disposable cannula comes sealed in a pouch and is Ethylene Oxide (ETO) sterilized (Figure 2). ESPX5 contains an individually wrapped cannula and an individually pouched endometrial suction pipet curette; ESPY5 contains only the individually wrapped cannula.

Check that the expiration date has not passed before use.


Figure 2: Disposable Cannula in pouch

- Packing Materials: Save the packing materials for potential future transportation and storage of the Endosee System®.
- Inspection: Inspect all components for damage during shipment, or discrepancies upon arrival. Ensure there are no unintended rough surfaces, sharp edges or protusions which may cause harm.

4.2 Power Requirements

The Endosee® System operates on DC voltage, from an internal battery source. The Display Module must be fully charged before the initial use.

- Connect the Display Module to the Docking Station, as shown below:

Figure 3: Display Module connected to Docking Station.

- Attach power supply cord to the Docking Station.
- Plug the power supply to a wall outlet.
- Continue to charge until the LED indicator on the Docking Station displays a steady green light, indicating the Display Module is fully charged. Once it is charged, the Display Module is ready for use.

4.3 Environmental Requirements

All of the environmental requirements listed in the Warnings and Precautions section must be met in the location where the Endosee® System is intended to be operated for patient use. If ANY of the seen environmental parameters are present, the Endosee® System should not be operated until modifications to the environmental conditions are met.

5 Description of Components

5.1 Product Description

The Endosee® System is comprised of the following components:

5.2 Reusable Handheld Display Module

The Display Module contains electronics including a power on/off button, a video processor, a display unit (LCD display), a rechargeable battery, battery management electronics, microcontrollers, storage for images and video, and firmware. The device uses a rechargeable battery. The battery is charged when the Display Module is connected to the Docking Station and the Docking Station is plugged into a wall outlet. A power supply is provided with the device. Figure 4 shows the reusable Display Module. The detailed battery requirements are listed in section 5.5.

1. LCD Touch Screen Monitor
2. Power On/Off Button

Figure 4: Reusable Display Module on Docking Station

5.3 Docking Station

The Docking Station has dual functions. It allows the battery to be charged using the power supply and allows video recordings and single frame photographs to be transferred via USB Cable from the internal memory of the Display Module to an external computer. Recordings and single frame photographs can also be downloaded directly from the Display Module via Wi-Fi.

WARNING: Use only the power supply provided with the device. Use of other chargers may damage the device and will void the warranty.

WARNING: Charge the battery at room temperature.


WARNING: Do not charge the battery near a heat source.

Figure 5: The Display Module on the Docking Station

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The disposable Cannula contains a miniature CMOS camera and light-emitting diode (LED) illumination module at the tip. The cannula connects to the Handheld Display Module through an electrical connector for: DC power, image data transfer and removal and control of the CMOS camera and illumination module. The cannula has an inner channel for the infusion of fluid. Fluid infusion is achieved by attaching either a pre-filled syringe or a pressurized bag of distending medium to the inflow port through an extension tube or a tubing set. The outflow port is at the distal tip of the cannula and empties out around the tip to infuse fluid into the uterine cavity or female bladder. The cannula has a pre-curved tip section at an angle of $20^{\circ} \pm 5^{\circ}$. The cannula tip section should not be altered. The cannula maximum temperature will not exceed 41°C . Figure 6 shows the cannula.

1. Handheld Display Module Connector
2. Inflow Port
3. Tip Section
4. CMOS Camera & LED Illumination Module

Figure 6: Disposable Cannula

5.5 Device Classification, Technical and Safety Specification

5.5.1. Technical Specifications

DISPLAY MODULE

Diagonal Size of Display Area:	125 mm (4.9 in)
Thickness:	22 mm (0.9 in)
Horizontal:	136 mm (5.4 in)
Height:	91 mm (3.6 in)
Weight (cannula, display):	281 g (0.6 lbs)

STERILE, SINGLE-USE CANNULA

Cannula Working Length:	278 mm (11 in)
Tip Largest Outside Dimension:	4.3 mm (0.17 in)
Pre-fixed tip angle:	20° ± 5°

POWER ADAPTER REQUIREMENTS

Input Power Range:	100-240V AC/0.2-0.4A
Frequency:	50/60Hz
Output:	5.99V DC/2.0A
Cable Length:	4 ft.
Charging Connector:	AC Adapter incorporated into Docking Station

BATTERY SPECIFICATIONS

Battery Type:	Lithium-ion rechargeable battery
Capacity:	3.7V rating, 1320mAh
Maximum Voltage when Battery is fully charged:	4.2 V
Minimum working voltage is:	≤ 3.0V
Full Charge Time:	Approx. 2 hours
Operating Time:	1 hour

SITE REQUIREMENTS — OPERATING CONDITIONS

Operating Temperature:	Environmentally Controlled Room Temperature
-------------------------------	---

6 Basics of Operation and Procedure

6.1 Cautions for Handling and Use

The following cautions should be strictly observed in the handling and use of the Endosee System.

WARNING: *The Endosee® System is for use only by physicians with adequate training in hysteroscopy/cystoscopy.*

6.2 Fluid Delivery System

Fluid infusion is achieved by attaching either a pre-filled syringe or a pressurized bag of distending medium to the inflow port of the cannula through an extension tube or a tubing set.

WARNING: *When a liquid medium is used, strict fluid intake and output surveillance should be maintained. Intrauterine instillation exceeding 1 liter should be followed with great care to reduce the possibility of fluid overload.*

WARNING: *Potential complications of Continuous Flow Hysteroscopy include:*

- Hyponatremia
- Hypothermia
- Uterine perforation resulting in possible injury to bowel, bladder, major blood vessels and ureter
- Pulmonary edema
- Cerebral edema
- Infection
- Bleeding
- Pain

The use of normal saline as a distending medium, and limiting the infused volume to less than 1000 mL is recommended to decrease the risk of the above complications. Intrauterine distension can usually be accomplished with pressures in the range of 35-75 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary

CAUTION: Do not allow any fluids or substances to be spilled in or around the cannula connector on the Handheld Display Module, as damage to the electronics may occur.

6.3 Cannula attachment and removal to the Display Module

In preparation for use, attach the cannula to the electrical connector on the Display Module per the following:

1. Peel apart the proximal end of the package and remove the cannula.
2. Open Cannula Handle to the “Down” position as shown as shown in Figure 7. Note that the Display Module will not connect to the Cannula unless the handle is in the “Down” position.
3. Align the connector faces of the Display Module and Cannula, and slide the Display Module onto the Cannula until the connectors are fully attached together.

CAUTION: Do not misalign the connectors between the cannula and the Handheld Display Module during assembly, as damage may occur.


Figure 7: Cannula and Display Module Attachment

6.4 Power Button Function

The function of the Display Module Power Button is described below.

WARNING: *Inspect the integrity of the unit and condition before powering-on the fully assembled unit. Do not use the system if visual signs of external damage are observed.*

Figure 8: Power Button

Touch Function	Action
Power On	Press-and-Release; Endosee splash screen will appear
Power Off	Press-and-Hold until audible beeps occur (approximately 1 second); Endosee splash screen will appear

6.5 Video and Picture (Single Frame) Button

The Video/Single Frame Buttons, located on both sides of the cannula, allow user to capture/save still images and video images. The image capture operations are described below.


Figure 9: Pictures (Single Frame) or Video Button


Video/Single frame Touch Function	Action
Capture Still Image	Press-and-Release
Capture Video Clip	Press-and-Hold: 1 second
Capture Still Image while Video is Recording	Press-and-Release while video is recording
Stop Video Capture	While system is capturing video, press button one more time

6.6 Touch Screen Functions

Touch Screen Description/Touch Screen Action


HOME MENU

Home menu has a battery power usage indicator and (4) four soft buttons each allowing a discrete function. Pressing each soft button enters a sub menu for a specific function. The instructions for touch screen user menus are described in this section.



Battery Power Usage Indicator

Battery usage indicator provides the battery power usage status. When battery power usage reaches 30%, the system will prompt user to charge the Display Module. When battery power usage reaches 20%, the system will automatically shut down.






Start Menu

- Upon entering Start Menu, the user has a choice of entering Exam mode for a New Patient or Saved Patient.
- If New Patient, press the desired numbers and/or letters to enter Patient ID. The numbers are displayed on the first screen. If letters are desired, select the soft button labeled “ABC” on the left side of the screen. User may enter up to twelve (12) characters.
- If Patient ID is not desired, a press “Back” to return to Exam Screen, and no Patient ID is entered.
- If Patient ID is entered incorrectly, press to cancel previously entered number.



- Press “Next” to confirm Patient ID and start the procedure. The newly created Patient ID will be displayed on the left side of the Exam screen.
- Verify the Display Module screen displays a live image in the correct orientation.
- To record images or video, follow instructions laid out on Section 6.5. The left side of the Exam screen will display a thumbnail of the previous image/video recorded.
- To exit Exam mode, touch “End Exam” on the bottom left of the Exam screen. The user will be returned to the Home menu.
- If “Saved Patient” is selected, the user will be taken to a screen with previously entered patient IDs and associated dates. To scroll through this list, touch the up or down arrows.
- Once the desired patient is selected, touch the “Next” button to enter Exam mode for that Patient. The selected patient ID will be displayed on the left side of the Exam screen.
- To record images or video, follow instructions laid out in Section 6.5. The left side of the Exam screen will display a thumbnail of the previous image/video recorded.
- To exit Exam mode, touch “End Exam” on the bottom left of the Exam screen. The user will be returned to the **Home menu**.

Touch Screen Description/Touch Screen Action



Preview Menu

- Pressing the **Preview** soft button in the **Home Menu** starts to display a real-time image. A pop-up message notifying the user of the file location for captured image/video will appear.
- Verify the Display Module screen displays a live image in the correct orientation
- Press one or both of the **Video/Single Frame Buttons** on the sides of the handheld Display Module to take a picture or record a video. Refer to 6.5 for detailed instructions.
- Press the “End Exam” soft button to return to the **Home menu**.

Playback Menu


- From the Home menu, press the “Playback” soft button on the display.
- Use the yellow up and down arrows to select the procedure of interest. The procedure can be found by the date or Patient ID that was entered.
- Once the procedure of interest is highlighted, press the “Playback” button in the bottom right corner of the display to review the procedure images and video.
- In the playback screen, the patient ID will be displayed in the bar on the left side of the screen.
- Upon first entering the Playback screen, the first image or video taken during the procedure will be displayed.
- Use the right and left arrows to scroll through the images and videos taken during the procedure.
- To play a video, press the “|▶” button at the bottom of the screen. Press it again to pause/restart the video.
- Press the “back” soft button on the bottom left of the screen to return to the “Playback” screen.
- Press the “back” soft button on the bottom left of the screen to return to the **Home menu**

Touch Screen Description/Touch Screen Action

Settings Menu

- From the Home menu, press the “Settings” soft button on the display. There are five menus in the Settings Menu: Date/ Time, Time Zone, WiFi, Pass Phrase and Service.

Time Zone


The time zone must be set by the user according to their locale. This is required as a network-connected device. These are the time zones available:

Eastern (USA)	Alaska (USA)	Hawaii (USA)	Eastern Standard (Canada)
Pacific (USA)	Atlantic Standard (Canada)	Eastern (Canada)	Mountain Standard (Canada)
Atlantic (Canada)	Central Standard (Canada)	Mountain (Canada)	Central (Mexico)
Central (Canada)	Central Standard (Canada)	Central (Mexico)	Guatemala
Pacific (Canada)	Eastern Standard (Mexico)	Guatemala	Mountain (Mexico)
Mountain Standard (Mexico)	Pacific (Mexico)	Belize	Honduras
El Salvador	Pacific (Mexico)	Mountain Standard (USA, AZ)	Costa Rica
Panama	Nicaragua	Newfoundland (Canada)	
Central (USA)	Mountain (USA)		

Date and Time Setup

From the **Home** menu, press the “Settings” soft button on the display. There are four menus in the Settings Menu: Date/Time, WiFi Mode, Wifi Passphrase, and Service.

- Use the yellow up and down arrows to highlight Date/Time. Once highlighted, press the “Next” button in the bottom right corner of the display to go to Date/Time screen.
- In the Date/Time screen, use the up and down arrows to set the date (MM/DD/YYYY) and time (HH:MM).
- Press “AM” or “PM” to complete the time of day setting.
- Once the date and time have been confirmed as accurate, press the “Save” soft button in the bottom right corner of the screen to return to the Settings screen.
- Press the “cancel” soft button in the bottom left of the screen if the Date/Time screen has been entered unintentionally.
- Press the “Back” soft button on the bottom left of the screen to return to **Home menu**.

WiFi Passphrase

Endosee features wireless access to retrieve videos and images. As in all wireless systems, proper care must be taken to insure the privacy and security of patient data. Wireless access to your Endosee is disabled by default. If you chose to enable it, the images and videos stored on the device are protected by a WIFI Passphrase or password of your choosing.

The first time you enable wireless access you must choose a password. This password will be retained unless you choose to change it from the menu.



The Endosee can store images and videos of exams. If you enter a name or ID for the patient, this field will also be stored with the image or video. This type of information is considered protected or private, so choose your password with care, and make sure that it is not generally available. Anyone who knows the password and has proximity to the Endosee can retrieve the exam data. It is good practice to change the password if an employee with access leaves the practice.

Password guidelines:

- A password should be 8 characters long minimum and can consist of a combination of numbers and letters.
- Longer passwords are generally better.
- The password should not be easily guessable. Random collections of characters are better than common words or obvious numeric sequences. Here are a few examples of poor passwords that are easily guessed:
 - Passwords based on dates
 - Passwords based on addresses
 - Passwords that form a pattern on a keypad
 - Passwords incorporating names or nicknames
 - Passwords incorporating personal data or preferences
- Don't reuse the same password for multiple devices or systems.
- Safeguard your password by storing it in a dedicated Password keeper.

Touch Screen Description/Touch Screen Action

- In the **Settings menu**, use the yellow up and down arrows to highlight WiFi Passphrase. Once highlighted, press the “next” button in the bottom right corner of the display to go to **WiFi Passphrase Screen**.
- Enter a WiFi Passphrase (note: a minimum of 8 characters is required). To create/set a WiFi Passphrase. Once entered, press the “next” button in the bottom right corner of the display and re-enter passphrase when prompted
- Press the “Next” button to save the passphrase and return to the **Settings Menu**. Note this passphrase can be reset by the user by repeating the previous step.

WiFi Mode

- In the Settings Menu use the yellow up and down arrows to highlight **WiFi Mode**. Once highlighted, press the “next” button in the bottom right corner of the display to go to **WiFi Mode Screen**.
- To enable WiFi, set **WiFi Mode** to “Temporary On” and the WiFi will remain on until the Display Module unit restarts. Or set Wifi Mode to “Always On” and the WiFi will remain on. Press the “Save” soft button to save your Wi-Fi Mode and press “Back” to return to the **Settings Menu**.
- To disable WiFi, set WiFi Mode to “Off” and the WiFi will remain off. Press the “Save” soft button to save your **Wi-Fi Mode** and press “Back” to return to the **Settings Menu**.

Service Menu

- The service screen is for use by CooperSurgical staff only.

6.7 Patient Examination Procedure

The following steps describe the process for using the Endosee® System Hysteroscope.

CAUTION: *A fully charged battery will provide at least 1 hour of continuous operating time under normal operating conditions. The Display Module should be fully charged before initial use. It is recommended to fully charge the Display Module after each procedure.*

CAUTION: *Refer to battery power indicator on the touch screen Display Module to determine remaining battery power.*

- Ensure that the Display Module is thoroughly cleaned and disinfected using the procedures, as described in Section 9.
- Ensure that the Display Module is charged, as described in Section 6.10.
- Remove the sterile cannula from the sterilized pouch.
- Inspect the cannula as described in Section 4.1.
- Attach the cannula to the Display Module using the procedures, as described in Section 6.3.
- Attach the fluid tubing to the cannula, as described in Section 6.2.
- Press and Release POWER ON/OFF button.
- The light source at the distal tip of the cannula will now be illuminated. If the tip of the cannula does not illuminate, system hardware troubleshooting procedures should be performed, as described in Section 12.
- Flush the fluid through the cannula until all air bubbles have cleared from the tubing.
- Using the Display Module touch screen, enter the Patient ID via the “Start” menu, or touch the Preview soft button as described in Section 6.6.

- Verify the Display Module Screen displays a live image (rather than a stored one) in the correct orientation.
- The Endosee® System is now ready for use.
- The cannula comes with a pre-curved tip at an angle of $20^{\circ} \pm 5^{\circ}$. The curve of the tip is similar to that of most conventional cervical dilators.
- The cannula is inserted through the cervical canal or uretha, with fluid flowing, while viewing the image, as is standard procedure for introduction of a hysteroscope or cystoscope.

6.8 Removal of the System

- After the procedure is completed, remove the cannula, and press and hold the Display Module POWER ON/OFF.
- Hold the Display Module such that the cannula tip is pointed towards the floor. Gently pull the cannula straight out from the Display Module. The cannula will detach from the Display Module.

CAUTION: When detaching the cannula from the Display Module, be sure that fluid does not contact the connector in the Display Module.

6.9 Disposal of Biohazardous Material

WARNING: The used cannula and/or disposable Separation Cable must be disposed of as biohazardous waste according to the safety guidelines of user facility/institution.

6.10 Battery Charging

The Endosee® System has an internal battery that powers the unit during use. A fully charged battery will provide at least 1 hour of continuous operating time under normal operating conditions.

To recharge the battery, follow the procedure in Section 4.2

7 External Interface Information

Docking Station

The Docking Station provides the capability for battery charging and image uploading.

Figure 10: Display Module on the Docking Station

DOWNLOADING STORED IMAGES

The unit should be kept in a secure location in order to ensure patient privacy.

Download to a PC via Docking Station

1. Plug Charging/Docking Station charge cable to Display Module.
2. Connect the USB cable to the computer and Charging/Docking Station.
3. Turn on Display Module power
4. Boot up computer.
5. Click Start, and then click My Computer.
6. The "ENDOSEE" device shows up as a removable disk in "My Computer." Open the ENDOSEE device to find a folder named with a patient ID.
7. Click folder named with patient ID.
8. The files in the exam date folder are recorded pictures and names end in "jpg."
9. Click folder named Preview. The files in the Preview folder are named by the date and time and ends in "jpg."
10. Copy those files by selecting them and copying them. The keyboard shortcut (control-A) marks all the photos in a folder and the keyboard shortcut (control-C) copies them to the computer's clipboard.

11. Navigate to the folder where you want to save patient pictures and paste (control-V) your photos into place.
12. After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
13. Do not erase photos from your memory card until you have confirmed that you saved them on the hard drive. It is recommended to back up file to another device.

Wireless Technology and Communication:

- The Endosee supports wireless standards of IEEE 802.11 a/b/g/n with a range of at least 10 feet. It can only connect to devices supporting WPA2-PSK AES protocol. It uses this to establish its own encrypted network and cannot connect to a local network or the internet. Because it cannot connect to external networks and does not contain removable media, the only cybersecurity responsibility imposed on the end user is limited to ensuring proper Endosee network password handling.
- The Endosee wireless access feature is only to retrieve patient videos and images after the procedure is complete. During the patient procedure, the Endosee device cannot be controlled or configured, or operated by wireless access.
- The Endosee device also provides USB connection to a computer for retrieving patient videos and images after the procedure is complete, which is a redundant function of the wireless access feature, in case the wireless access has a failure.

Download to a PC via WiFi

1. Turn on Display Module power
2. Ensure WiFi setting is Enabled (see section 6.6 – Touch Screen Functions – WiFi Passphrase)
3. Boot up computer
4. Check available networks on the computer
5. Select the Endosee Device. Endosee Device should be identified as “ENDOSEE” followed by a serial number for the device unit.
6. Type in the passphrase previously entered in section 6.6 when prompted
7. ENDOSEE Device should be connected wirelessly to the computer

8. Open Windows File Explorer. Depending on your PC's network configuration, an Endosee Device may show up under "Network" in the File Explorer window. If no ENDOSEE device is detected, enter \\192.168.0.1 to the file explorer address bar.
9. Open the ENDOSEE device to find a folder with named with a Patient ID.
10. The files in the exam date folder are recorded pictures and names end in "jpg."
11. Click folder named Preview
12. The files in the Preview folder are named by the date and time and ends in "jpg."
13. Copy those files by selecting them and copying them. The keyboard shortcut (command-A) marks all the photos in a folder and the keyboard short cut (command-C) copies them to the computer's clipboard. Navigate to the folder where you want to save patient pictures and paste (command-V) your photos into place.
14. After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
15. Do not erase photos from your memory card until you have confirmed that you saved them on the hard drive.

Download to a Mac via Docking Station

1. Plug Charging/Docking Station charge cable to Display Module.
2. Connect the USB cable to the computer and Charging/Docking Station.
3. Turn on Display Module power.
4. Boot up computer.
5. Click the Finder on the application dock.
6. The "ENDOSEE" device shows up as a removable device. Open the ENDOSEE device to find a folder named with a Patient ID.
7. The files in the exam date folder are recorded pictures and names end in "jpg."
8. Click folder named Preview

9. The files in the Preview folder are named by the date and time and ends in “jpg.”
10. Copy those files by selecting them and copying them. The keyboard shortcut (command-A) marks all the photos in a folder and the keyboard shortcut (command-C) copies them to the computer’s clipboard.
11. Navigate to the folder where you want to save patient pictures and paste (command-V) your photos into place.
12. After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
13. Do not erase photos from your memory card until you have confirmed that you saved them on the hard drive.

Download to a Mac via WiFi

1. Turn on Display Module power
2. Ensure WiFi setting is Enabled (see section 6.6 – Touch Screen Functions – WiFi Passphrase)
3. Boot up computer
4. Check available networks on the computer
5. Select the Endosee Device. Endosee Device should be identified as “ENDOSEE” followed by a serial number for the device unit.
6. Type in the passphrase previously entered in section 6.6 when prompted
7. ENDOSEE device should be connected wirelessly to the computer
8. Click on Finder on the application dock and locate the Endosee Device under “Shared” devices. If no ENDOSEE device is detected, enter \\192.168.0.1 to the finder address bar.
9. Open the ENDOSEE device to find a folder named with a patient ID.
10. The files in the exam date folder are recorded pictures and names end in “jpg.”
11. Click folder named Preview

12. The files in the Preview folder are named by the date and time and ends in “jpg.”
13. Copy those files by selecting them and copying them. The keyboard shortcut (command-A) marks all the photos in a folder and the keyboard shortcut (command-C) copies them to the computer’s clipboard.
14. Navigate to the folder where you want to save patient pictures and paste (command-V) your photos into place.
15. After you copy files or folders to the hard drive, it is useful to view the hard drive to confirm that the files are copied.
16. Do not erase photos from your memory card until you have confirmed that you saved them on the hard drive.

** It is recommended to save the last month of cases in two storage locations, as a back up.*

8 EMC Compliance Information

8.1 Guidance and Manufacturer’s Declaration – Electromagnetic emissions

The Endosee System is intended for use in the electromagnetic environment specified below. The customer or the user of the Endosee System should assure that it is used in such an environment.		
Professional Healthcare Facility Environment		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Endosee System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

8.2 Guidance and Manufacturer's Declaration – Electromagnetic immunity

<p>The Endosee System is intended for use in the electromagnetic environment specified below. The customer or the user of the Endosee System should assure that it is used in such an environment.</p>		
<p>Professional Healthcare Facility Environment</p>		
<p>Immunity Test</p>	<p>IEC60601 Test and Compliance</p>	<p>Electromagnetic environment – guidance</p>
<p>Electrostatic discharge (ESD) IEC 61000-4-2</p>	<p>+/- 8 kV contact +/- 15 kV air</p>	
<p>Radiated RF EM Fields IEC 61000-4-3</p>	<p>Professional Health Care 3V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz</p>	
<p>Proximity fields from RF wireless communications equipment IEC 61000-4-3</p>	<p>9 V/m 710,745,780 MHz Pulse modulation 217 Hz</p> <p>9 V/m 5240,5500, 5785 MHz Pulse modulation 217 Hz</p> <p>27 V/m 385 MHz Pulse modulation 18Hz</p> <p>28 V/m 450 MHz</p>	<p>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Endosee® System, including cables specified by the manufacturer.</p>
	<p>FM +/- 5 kHz deviation; 1 kHz sine</p> <p>28 V/m 810, 870, 930 MHz Pulse modulation 18 Hz</p> <p>28 V/m 1720,1845,1970, 2450 MHz Pulse modulation 217 Hz</p>	<p>Otherwise, degradation of the performance of this equipment could result.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Professional Healthcare Facility Environment		
Immunity Test	IEC60601 Test and Compliance	Electromagnetic environment – guidance
Electrical Fast Transient/Burst IEC 61000-4-4	+/- 2 kV for ac power supply lines +/- 1 kV for input/output lines 100 kHz repetition frequency	
Surge IEC 61000-4-5	+/- 0.5, 1 kV line(s) to line(s) +/- 0.5, 1, 2 kV line(s) to earth	
Conducted disturbances induced by RF fields IEC 61000-4-6	Professional Health Care 3V 0.15-80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1kHz	
Voltage dips IEC 61000-4-11	0 % U_T for 0.5 cycle 0, 45, 90, 135, 180, 225, 270 and 315° 0 % U_T for 1 cycle 0° 70 % U_T for 25/30 cycles 0°	
Voltage Interruptions IEC 61000-4-11	0 % U_T for 250/300 cycles 0°	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 or 60 Hz	
NOTE U_T is the a.c. mains voltage prior to application of the test level.		

9 Maintenance

9.1 Cleaning / Disinfecting Instructions

After each patient use, follow the combined cleaning/ disinfecting instructions provided below for the Display Module.

WARNING: Follow all cleaning procedures provided.

CAUTION: No portion of the Display Module should be opened for cleaning or disinfecting.

9.2 Intermediate Level Combined Cleaning/Disinfecting Instructions

- After each patient use, power off and unplug all electric connections.
- Place Display Module on a clean, dry, disinfected or sterile basin or surface.

CAUTION: Throughout the cleaning/ disinfection process, orient the connector in a manner that prevents liquid from entering the electrical connection area of the connector (e.g., with the opening to the electrical connections facing upward). Do not wet the opening or inside of the connector housing the electrical connection components.

CAUTION: When cleaning/disinfecting the LCD screen of the Display Module, pay careful attention to the seams around the edge of the display and the back of the unit, including seam around Power ON/OFF Button.

- Follow the instructions below to clean and disinfect the Display Module. These instructions were validated using CaviWipes® sanitizing wipes. If you choose a wipe not listed here, you must validate it before use.
- Using a fresh, moist disinfectant wipe, thoroughly wipe the front and back of the LCD screen of the Display Module. Discard the wipe. Repeat with a fresh wipe for a minimum of thirty (30) seconds. Discard the wipe.
- Using a fresh, moist wipe, thoroughly wipe the grip section of the Display Module. Discard the wipe. Repeat with a fresh wipe for a minimum of ten (10) seconds. Discard the wipe.
- Using a fresh, moist wipe, thoroughly wipe the seams of the Display Module for one (1) minute.
- Using a fresh, moist wipe, thoroughly wipe all surfaces of the Display Module for a minimum of ten (10) seconds. Discard the wipe.
- Allow the wiped and wetted Display Module to stand for a minimum of three (3) minutes.
- At the end of the three (3) minutes, use a sterile, lint-free wipe saturated with sterile 70% IPA to thoroughly wipe the entire Handheld Monitor for a minimum of thirty (30) seconds.
- Allow the device to air dry. Once dry, visually inspect the Display Module for any residual soil. Repeat cleaning/disinfection instructions if visible soil is present.
- Place the cleaned/disinfected device in an appropriate dry storage area.

9.3 Disposable Cannula Handling/Care

The cannula is a single-use, disposable item. As such, no formal cleaning or disinfection procedures apply to the cannula.

9.4 Maintenance

The Endosee® System has no user - serviceable parts. Note: The Endosee® System requires no calibration or service in the field. All electronics are tested and validated for performance at the manufacturer. If the performance of the Endosee Cannula is in doubt, contact CooperSurgical, Inc. Customer Support for assistance.

Maintenance Intervals

Basic preventive maintenance of the Endosee® System is an important function that should be performed on schedule to ensure safe and effective operation. If the integrity or performance of any portion of the Endosee® System is in doubt, the system should NOT be used until the issue is resolved.

The following maintenance intervals are recommended:

Maintenance Type	Interval	Personnel
General Inspection	Before each Patient Use	Users
Safety Label Integrity	As Required	All

10 Storage & Shipping

The Endosee® System is shipped in a protective box. Do not dispose of protective box or any other shipping materials. These items should be retained for future storage or transport of the equipment.

CAUTION: Do not attempt to transport or ship the Endosee® System without using proper packaging to protect the product. Store at room temperature.

11 Technical Assistance

For technical information, assistance or to order parts or manuals, contact CooperSurgical, Inc. Customer Support Department, or your local CooperSurgical, Inc. representative.

Information regarding system serial number is available from Customer Support.

95 Corporate Drive • Trumbull, CT 06611 • USA
Phone: (800) 243-2974 • Fax: (800) 262-0105
www.coopersurgical.com

12 Troubleshooting Guide

The information in this section is intended to provide simple steps that can be performed by a user for identifying the primary cause, and possible simple solutions that can be resolved on-site, to basic problems that may be encountered while operating the Endosee® System.

Any issues determined to be beyond the scope of the basic user troubleshooting steps provided in this User Manual should be communicated to the CooperSurgical, Inc. Customer Support department.

WARNING: *The Endosee® System contains no user-serviceable components within the enclosure. To avoid electrical shock or damage, the unit must not be disassembled. Doing so will void the warranty.*

WARNING: *Do not disassemble or alter the battery in any way.*

WARNING: *Battery is not user replaceable.*

Problem	Test	Action
(1) When switched ON, there is no display on screen	Is the cannula tip LED illuminated?	<ul style="list-style-type: none"> • If NO, first ensure cannula is fully inserted in the Display Module. If still no illumination/ screen display, switch to a new cannula. • If YES, contact CooperSurgical, Inc. Customer Service.
(2) Poor picture quality	N/A	<ul style="list-style-type: none"> • Clean cannula tip with a sterile, clean, soft wipe. • Ensure cannula is fully inserted in the Display Module.
(3) Image on screen flickers or has lines across it	N/A	<ul style="list-style-type: none"> • Ensure cannula is fully inserted in the Display Module.
(4) Loose component, poorly fitting connection	N/A	<ul style="list-style-type: none"> • Do not use device.
(5) <i>Image displays on half of the screen</i>	N/A	<ul style="list-style-type: none"> • Disconnect and reconnect the cannula.
(6) Touchscreen is unresponsive	N/A	<ul style="list-style-type: none"> • Press and hold the Power Button for 5 seconds and restart the unit.
(7) No Data Transfer/File access over USB	N/A	<ul style="list-style-type: none"> • Instruct hospital IT department to configure the computer to access files via USB 2.0.
(8) No Data Transfer/File access over WiFi	N/A	<ul style="list-style-type: none"> • Instruct hospital IT department to configure the computer to access files over WiFi

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© 2019 CooperSurgical, Inc.
CaviWipes® is a registered trademark of Metrex Research Corporation.
CooperSurgical, Inc. • 95 Corporate Drive • Trumbull, CT 06611 USA